


## Product Datasheet

### RGS17 Antibody (orb1254102)

<b>Description</b>	RGS17 Antibody
<b>Species/Host</b>	Gallus
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	41-100
<b>Target</b>	RGS17
<b>Preservatives</b>	Phosphate-Buffered Saline. No preservatives added.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	RGS17 antibody can be stored at 4°C for short term (weeks). Long term storage should be at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
<b>Note</b>	For research use only
<b>Application notes</b>	RGS17 antibody can be used for the detection of RGS17 by Western blot, may also work for IHC and ICC.
<b>Clonality</b>	Polyclonal
<b>MW</b>	24.4 kDa (calculated)
<b>Uniprot ID</b>	<b>Q9UGC6</b>
<b>NCBI</b>	<b>NP_036551.3</b>
<b>Dilution Range</b>	RGS17 antibody can be used for the detection of RGS17 by Western blot, may also work for IHC and ICC.
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot:  
orb1254102  
dilution: 1:2,  
...