



Product Datasheet

COL10A1 Antibody (orb1254362)


Description	COL10A1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 521-680 of human COL10A1 (NP_000484.2).
Target	COL10A1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IHC: 1:100 - 1:200
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q03692
Dilution Range	IHC: 1:100 - 1:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...