


Product Datasheet

P2RY11 Antibody (orb1254649)

Description	P2RY11 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant protein of human P2RY11
Target	P2RY11
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IHC: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96G91
Dilution Range	IHC: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...