


Product Datasheet

MT2A Antibody (orb1254713)

Description	MT2A Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant protein of human MT2A.
Target	MT2A
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02795
Dilution Range	IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...