



Product Datasheet

SCGB1A1 Antibody (orb1254728)


| | |
|----------------------------|---|
| Description | SCGB1A1 Antibody |
| Species/Host | Rabbit |
| Reactivity | Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 22-91 of human SCGB1A1 (NP_003348.1). |
| Target | SCGB1A1 |
| Preservatives | Liquid |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Store at -20˚°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Application notes | IHC: 1:50 - 1:200IF: 1:50 - 1:200 |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P11684 |
| Dilution Range | IHC: 1:50 - 1:200IF: 1:50 - 1:200 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...



Immunofluorescence analysis of mouse lun...