



Product Datasheet

KDM7A Antibody (orb1254817)

Description	KDM7A Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 417-735 of mouse KDM7A (NP_001028602.2).
Target	KDM7A
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HCT116 ce...



Immunofluorescence analysis of U2OS cell...