


## Product Datasheet

### FURIN Antibody (orb1255079)

<b>Description</b>	FURIN Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A synthesized peptide derived from human Furin
<b>Target</b>	FURIN
<b>Preservatives</b>	Liquid
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4&#730°C short term. Aliquot and store at -20&#730°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: 1:500 - 1:2000
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>MW</b>	Predicted: 87 kDa Observed: 90 kDa
<b>Uniprot ID</b>	<a href="#">P09958</a>
<b>NCBI</b>	<a href="#">P09958</a>
<b>Dilution Range</b>	WB: 1:500 - 1:2000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of extracts of var...