



Product Datasheet

CCR5 (phospho-Ser349) antibody (orb125523)

Description	Rabbit polyclonal antibody to CCR5 (phospho-Ser349).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P51681
Dilution Range	WB: 1:500-1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Western blot analysis of extracts from A...