



Product Datasheet

Mt-Cyb Antibody (orb1255582)

Description	Mt-Cyb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse MT-CYB (NP_904340.1).
Target	Mt-Cyb
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:200 - 1:1000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 43kDa
Dilution Range	WB: 1:200 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Western blot analysis of extracts of var...