


Product Datasheet

EBNA1BP2 Antibody (orb1255654)

Description	EBNA1BP2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of human EBNA1BP2 (NP_006815.2).
Target	EBNA1BP2
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:2000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 40kDa
Uniprot ID	Q99848
Dilution Range	WB: 1:1000 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...