



Product Datasheet

SART1 Antibody (orb1255762)

Description	SART1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human SART1 (NP_005137.1).
Target	SART1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IP: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 120kDa
Uniprot ID	O43290
Dilution Range	WB: 1:500 - 1:2000IP: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunoprecipitation analysis of 150 ug e...