



Product Datasheet

CHST13 antibody (orb125580)

Description	Rabbit Polyclonal to CHST13.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NET6
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from K...



Immunohistochemical analysis of formalin...