



Product Datasheet

PPP4R1 Antibody (orb1255885)


Description	PPP4R1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 324-647 of human PPP4R1 (NP_001035847.1).
Target	PPP4R1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IP: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 130kDa
Uniprot ID	Q8TF05
Dilution Range	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IP: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...