



Product Datasheet

ZP3 Antibody (orb1255891)


| | |
|----------------------------|--|
| Description | ZP3 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, WB |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 170-390 of human ZP3 (NP_001103824.1). |
| Target | ZP3 |
| Preservatives | Liquid |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Store at -20˚°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Application notes | WB: 1:500 - 1:2000IF: 1:50 - 1:200 |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | Observed: 47kDa |
| Uniprot ID | P21754 |
| Dilution Range | WB: 1:500 - 1:2000IF: 1:50 - 1:200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of SKO...



Immunofluorescence analysis of rat oopho...



Immunofluorescence analysis of rat oopho...