



Product Datasheet

NSUN6 Antibody (orb1256530)

Description	NSUN6 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-469 of human NSUN6 (NP_872349.1).
Target	NSUN6
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IP: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 52kDa
Uniprot ID	Q8TEA1
Dilution Range	WB: 1:500 - 1:2000IP: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunoprecipitation analysis of 200 ug e...