



Product Datasheet

TEAD1 Antibody (orb1256786)


Description	TEAD1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 135-215 of human TEAD1 (NP_068780.2).
Target	TEAD1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:50 - 1:200ChIP: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 50kDa
Uniprot ID	P28347
Dilution Range	WB: 1:500 - 1:2000IF: 1:50 - 1:200ChIP: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of MCF7 cell...



Chromatin immunoprecipitation analysis o...