



Product Datasheet

IL17F Antibody (orb1256978)

| | |
|----------------------------|--|
| Description | IL17F Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IF, WB |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 31-163 of human IL17F (NP_443104.1). |
| Target | IL17F |
| Preservatives | Liquid |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Application notes | WB: 1:500 - 1:2000IF: 1:50 - 1:100 |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | Observed: 18kDa |
| Uniprot ID | Q96PD4 |
| Dilution Range | WB: 1:500 - 1:2000IF: 1:50 - 1:100 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of var...



Immunofluorescence analysis of MCF7 cell...