



Product Datasheet

SSRP1 Antibody (orb1257028)


Description	SSRP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SSRP1 (NP_003137.1).
Target	SSRP1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IHC: 1:50 - 1:200IF: 1:10 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 98kDa
Uniprot ID	Q08945
Dilution Range	WB: 1:500 - 1:2000IHC: 1:50 - 1:200IF: 1:10 - 1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...