



Product Datasheet

CRABP2 Antibody (orb1257194)


Description	CRABP2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-94 of human CRABP2 (NP_001869.1).
Target	CRABP2
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IHC: 1:50 - 1:200IF: 1:20 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 15kDa
Uniprot ID	P29373
Dilution Range	WB: 1:500 - 1:2000IHC: 1:50 - 1:200IF: 1:20 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunofluorescence analysis of U2OS cell...