



Product Datasheet

XRCC5 Antibody (orb1257285)


Description	XRCC5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 543-732 of human XRCC5 (NP_066964.1).
Target	XRCC5
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200 IF: 1:20 - 1:100 IP: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 86kDa
Uniprot ID	P13010
Dilution Range	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200 IF: 1:20 - 1:100 IP: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...