



Product Datasheet

VAPB Antibody (orb1257531)


Description	VAPB Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-222 of human VAPB (NP_004729.1).
Target	VAPB
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 35kDa/47kDa
Uniprot ID	O95292
Dilution Range	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...