



Product Datasheet

SERPINA10 Antibody (orb1257654)


Description	SERPINA10 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-444 of human SERPINA10 (NP_001094077.1).
Target	SERPINA10
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:20 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 51kDa
Uniprot ID	Q9UK55
Dilution Range	WB: 1:500 - 1:2000IF: 1:20 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Hep...



Immunofluorescence analysis of A549 cell...



Immunofluorescence analysis of PC-12 cel...