


Product Datasheet

UBL4A Antibody (orb1257736)

| | |
|----------------------------|---|
| Description | UBL4A Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 1-157 of human UBL4A (NP_055050.1). |
| Target | UBL4A |
| Preservatives | Liquid |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Store at -20˚°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Application notes | WB: 1:200 - 1:1000 |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | Observed: 14-16kDa |
| Uniprot ID | P11441 |
| Dilution Range | WB: 1:200 - 1:1000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of var...