


Product Datasheet

THBD Antibody (orb1257747)

Description	THBD Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-500 of human THBD (NP_000352.1).
Target	THBD
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:200 - 1:2000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 110kDa
Uniprot ID	P07204
Dilution Range	WB: 1:200 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...