


Product Datasheet

SFTP A1 Antibody (orb1257941)

Description	SFTP A1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human SFTP A1
Target	SFTP A1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 35kDa
Uniprot ID	Q8IWL2
Dilution Range	WB: 1:500 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...