



Product Datasheet

CLDN3 Antibody (orb1258002)


Description	CLDN3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein of human CLDN3
Target	CLDN3
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 23kDa
Uniprot ID	O15551
Dilution Range	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...