



Product Datasheet

EIF4G2 Antibody (orb1258018)


Description	EIF4G2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 750-850 of human EIF4G2 (NP_001409.3).
Target	EIF4G2
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000 IHC: 1:100 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 102kDa
Uniprot ID	P78344
Dilution Range	WB: 1:500 - 1:1000 IHC: 1:100 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...