


Product Datasheet

DRD4 Antibody (orb1258021)

Description	DRD4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human DRD4
Target	DRD4
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 38kDa
Uniprot ID	P21917
Dilution Range	WB: 1:500 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...