



Product Datasheet

mTOR Antibody (orb1258164)


Description	mTOR Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human mTOR (NP_004949.1).
Target	mTOR
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:50 - 1:200IP: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 289kDa
Uniprot ID	P42345
Dilution Range	WB: 1:500 - 1:2000IF: 1:50 - 1:200IP: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of PC12 cell...