



Product Datasheet

FMO3 Antibody (orb1258400)

Description	FMO3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 263-532 of human FMO3 (NP_001002294.1).
Target	FMO3
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:10 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 68kDa
Uniprot ID	P31513
Dilution Range	WB: 1:500 - 1:2000IF: 1:10 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...



Immunofluorescence analysis of A-549 cel...