



Product Datasheet

CYP27B1 Antibody (orb1258759)

Description	CYP27B1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-195 of human CYP27B1 (NP_000776.1).
Target	CYP27B1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 63kDa
Uniprot ID	O15528
Dilution Range	WB: 1:500 - 1:2000IF: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of HeLa cell...