



Product Datasheet

MKI67 Antibody (orb1258827)

Description	MKI67 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein of human Ki67
Target	MKI67
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 358kDa
Uniprot ID	P46013
Dilution Range	WB: 1:500 - 1:2000IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.




Western blot analysis of extracts of NIH...



Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of U2OS cell...