



## Product Datasheet

### HOXD9 Antibody (orb1258841)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | HOXD9 Antibody   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | WB   |
| <b>Immunogen</b>           | A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HOXD9 (NP_055028.3). |
| <b>Target</b>              | HOXD9  |
| <b>Preservatives</b>       | Liquid   |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Concentration</b>       | batch dependent  |
| <b>Storage</b>             | Store at -20&#730°C. Avoid freeze / thaw cycles.   |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | WB: 1:500 - 1:2000   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <b>P28356</b>  |
| <b>Dilution Range</b>      | WB: 1:500 - 1:2000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Western blot analysis of extracts of SH-...



Western blot analysis of extracts of var...