


Product Datasheet

PPAP2B Antibody (orb1259359)

Description	PPAP2B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PPAP2B (NP_003704.3).
Target	PPAP2B
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 38kDa
Uniprot ID	O14495
Dilution Range	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...