


Product Datasheet

MDC1 (phospho-Ser513) antibody (orb125952)

Description	Rabbit Polyclonal to MDC1 (phospho-Ser513).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q14676
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...