



Product Datasheet

ROBO1 Antibody (orb1259569)


| | |
|----------------------------|---|
| Description | ROBO1 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 1452-1651 of human ROBO1 (NP_002932.1). |
| Target | ROBO1 |
| Preservatives | Liquid |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Application notes | WB: 1:200 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200 |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | Observed: 270kDa |
| Uniprot ID | Q9Y6N7 |
| Dilution Range | WB: 1:200 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of var...



Western blot analysis of extracts of var...



Immunofluorescence analysis of U-2 OS ce...