



Product Datasheet

MRPL39 antibody (orb125982)

Description	Rabbit Polyclonal to MRPL39.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NYK5
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from M...



Immunohistochemical analysis of formalin...