


Product Datasheet

CCDC83 Antibody (orb1259939)

| | |
|----------------------------|---|
| Description | CCDC83 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CCDC83 (NP_001273088.1). |
| Target | CCDC83 |
| Preservatives | Liquid |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Store at -20˚°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Application notes | WB: 1:500 - 1:2000 |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | Observed: 49kDa |
| Uniprot ID | Q8IWF9 |
| Dilution Range | WB: 1:500 - 1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of var...