



Product Datasheet

p47 phox (phospho-Ser345) antibody (orb126026)


Description	Rabbit polyclonal antibody to p47 phox (phospho-Ser345)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P14598
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells lys...



Western blot analysis of extracts from H...



Immunohistochemical analysis of formalin...