



Product Datasheet

PNLIPRP2 Antibody (orb1260692)


| | |
|----------------------------|--|
| Description | PNLIPRP2 Antibody |
| Species/Host | Rabbit |
| Reactivity | Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, WB |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 301-470 of human PNLIPRP2 (NP_005387.2). |
| Target | PNLIPRP2 |
| Preservatives | Liquid |
| Form/Appearance | Liquid |
| Concentration | batch dependent |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Application notes | WB: 1:500 - 1:1000IF: 1:50 - 1:200 |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | Observed: 54kDa |
| Uniprot ID | P54317 |
| Dilution Range | WB: 1:500 - 1:1000IF: 1:50 - 1:200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of rat...



Immunofluorescence analysis of rat pancreas...



Immunofluorescence analysis of mouse pancreas...