


Product Datasheet

TCL1B Antibody (orb1260735)

Description	TCL1B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human TCL1B (NP_004909.1).
Target	TCL1B
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 15kDa
Uniprot ID	095988
Dilution Range	WB: 1:500 - 1:1000IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...