



Product Datasheet

PTGS2 Antibody (orb1260868)


Description	PTGS2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 375-604 of human Cox2 (NP_000954.1).
Target	PTGS2
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:4000 IHC: 1:50 - 1:100 IF: 1:50 - 1:100 IP: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 80kDa
Uniprot ID	P35354
Dilution Range	WB: 1:1000 - 1:4000 IHC: 1:50 - 1:100 IF: 1:50 - 1:100 IP: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...