



Product Datasheet

TOP1 Antibody (orb1260871)


Description	TOP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TOP1 (NP_003277.1).
Target	TOP1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 120kDa
Uniprot ID	P11387
Dilution Range	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunoprecipitation analysis of 200 ug e...