


Product Datasheet

GNG7 Antibody (orb1261725)

Description	GNG7 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human GNG7 (NP_443079.1).
Target	GNG7
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 12kDa
Uniprot ID	O60262
Dilution Range	WB: 1:500 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...