



Product Datasheet

HDAC4 Antibody (orb1261911)


Description	HDAC4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 451-650 of human HDAC4 (NP_006028.2).
Target	HDAC4
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20˚°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000IHC: 1:50 - 1:100IF: 1:100 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 130kDa
Uniprot ID	P56524
Dilution Range	WB: 1:500 - 1:1000IHC: 1:50 - 1:100IF: 1:100 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...