



Product Datasheet

SGMS2 Antibody (orb1262078)

Description	SGMS2 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis of SGMS2 Antibody ...
Immunogen	This SGMS2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 338-365 amino acids from the C-terminal region of human SGMS2.	
Target	SGMS2	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	Western blot analysis of SGMS2 Antibody ...
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	Western blot analysis of SGMS2 using rab...
MW	42 kDa	
Uniprot ID	Q8NHU3	
NCBI	Q8NHU3	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	