

Product Datasheet

ABCD2 Antibody (orb1262130)

Description

ABCD2 Antibody


Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This ABCD2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 269-298 amino acids from the Central region of human ABCD2.

Western blot analysis of ABCD2 Antibody ...

Target

ABCD2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

83 kDa

Uniprot ID
[Q9UBJ2](#)
NCBI
[Q9UBJ2](#)
Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.