



Product Datasheet

FZD4 Antibody (orb1262201)

Description	FZD4 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	Western blot analysis of FZD4 Antibody i...
Tested Applications	WB	
Immunogen	This FZD4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 495-522 amino acids from the C-terminal region of human FZD4.	
Target	FZD4	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	Western blot analysis of FZD4 Antibody i...
Form/Appearance	Liquid	
Concentration	batch dependent	Western blot analysis of FZD4 Antibody i...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	60 kDa	
Uniprot ID	Q9ULV1	
NCBI	Q9ULV1	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	