



Product Datasheet

COQ5 Antibody (orb1262268)

Description	COQ5 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	Western blot analysis of COQ5 Antibody i...
Tested Applications	IHC-P, WB	
Immunogen	This COQ5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 213-241 amino acids from the Central region of human COQ5.	
Target	COQ5	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	COQ5 Antibody IHC analysis in formalin f...
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	37 kDa	
Uniprot ID	Q5HYK3	
NCBI	Q5HYK3	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	