
Product Datasheet

EIF4A2 Antibody (orb1262327)

Description

EIF4A2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IHC-P, WB

Immunogen

This EIF4A2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 333-360 amino acids from the C-terminal region of human EIF4A2.

Target

EIF4A2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


46 kDa

Uniprot ID
[Q14240](#)
NCBI
[Q14240](#)
Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot at 1:4000 dilution + MCF-7 ...



Western blot analysis of EIF4A2 Antibody...



Formalin-fixed and paraffin-embedded hum...